

PLL_P Antibody, Biotin conjugated



PACO42621

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Appears to be involved in myelination. Could also participate in ion transport events as addition of plasmolipin to lipid bilayers induces the formation of ion channels, which are voltage-dependent and K(+)-selective.

Gene ID:

PLL_P

Uniprot

Q9Y342

Synonyms:

Plasmolipin (Plasma membrane proteolipid), PLL_P, PMLP TM4SF11

Immunogen:

Recombinant Human Plasmolipin protein (1-35AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A