

TIMM10B Antibody, FITC conjugated



PACO42632

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Component of the TIM22 complex, a complex that mediates the import and insertion of multi-pass transmembrane proteins into the mitochondrial inner membrane. The TIM22 complex forms a twin-pore translocase that uses the membrane potential as the external driving force. In the TIM22 complex, it may act as a docking point for the soluble 70 kDa complex that guides the target proteins in transit through the aqueous mitochondrial intermembrane space.

Gene ID:

TIMM10B

Uniprot

Q9Y5J6

Synonyms:

Mitochondrial import inner membrane translocase subunit Tim10 B (Fracture callus protein 1) (FxC1) (Mitochondrial import inner membrane translocase subunit Tim9 B) (TIMM10B) (Tim10b), TIMM10B, FXC1 TIM9B TIMM9B

Immunogen:

Recombinant Human Mitochondrial import inner membrane translocase subunit Tim10 B protein (1-103AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A