

# MORC2 Antibody, HRP conjugated



PACO42647

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

Exhibits a cytosolic function in lipogenesis, adipogenic differentiation, and lipid homeostasis by increasing the activity of ACLY, possibly preventing its dephosphorylation. May act as a transcriptional repressor. Down-regulates CA9 expression.

**Gene ID:**

MORC2

**Uniprot**

Q9Y6X9

**Synonyms:**

MORC family CW-type zinc finger protein 2 (Zinc finger CW-type coiled-coil domain protein 1), MORC2, KIAA0852 ZCWCC1

**Immunogen:**

Recombinant Human MORC family CW-type zinc finger protein 2 protein (2-300AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A