

MORC2 Antibody, FITC conjugated



PACO42648

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Exhibits a cytosolic function in lipogenesis, adipogenic differentiation, and lipid homeostasis by increasing the activity of ACLY, possibly preventing its dephosphorylation. May act as a transcriptional repressor. Down-regulates CA9 expression.

Gene ID:

MORC2

Uniprot

Q9Y6X9

Synonyms:

MORC family CW-type zinc finger protein 2 (Zinc finger CW-type coiled-coil domain protein 1), MORC2, KIAA0852 ZCWCC1

Immunogen:

Recombinant Human MORC family CW-type zinc finger protein 2 protein (2-300AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A