

NOP58 Antibody, HRP conjugated



PACO42659

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Required for 60S ribosomal subunit biogenesis.

Gene ID:

NOP58

Uniprot

Q9Y2X3

Synonyms:

Nucleolar protein 58 (Nucleolar protein 5), NOP58, NOL5 NOP5

Immunogen:

Recombinant Human Nucleolar protein 58 protein (1-529AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A