

# NOP58 Antibody, FITC conjugated



PACO42660

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

Required for 60S ribosomal subunit biogenesis.

**Gene ID:**

NOP58

**Uniprot**

Q9Y2X3

**Synonyms:**

Nucleolar protein 58 (Nucleolar protein 5), NOP58, NOL5 NOP5

**Immunogen:**

Recombinant Human Nucleolar protein 58 protein (1-529AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A