

# AAR2 Antibody



PACO42662

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, IHC

**Recommended dilutions:**

ELISA:1:2000-1:10000, IHC:1:20-1:200

**Protein Background:**

spliceosomal tri-snRNP complex assembly.

**Gene ID:**

AAR2

**Uniprot**

Q9Y312

**Synonyms:**

Protein AAR2 homolog (AAR2 splicing factor homolog), AAR2, C20orf4

**Immunogen:**


Recombinant Human Protein AAR2 homolog protein (2-250AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---



Immunohistochemistry of paraffin-embedded human tonsil tissue using PACO42662 at dilution of 1:100.