

PRELID3B Antibody



PACO42670

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

mitochondrial intermembrane space, phosphatidic acid, transporter activity, phospholipid transport.

Gene ID:

PRELID3B

Uniprot

Q9Y3B1

Synonyms:

PRELI domain containing protein 3B (Protein slowmo homolog 2), PRELID3B, C20orf45 SLMO2


Immunogen:

Recombinant Human PRELI domain containing protein 3B protein (1-194AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human thyroid tissue using PACO42670 at dilution of 1:100.