

RBM7 Antibody, Biotin conjugated



PACO42689

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Possible involved in germ cell RNA processing and meiosis.

Gene ID:

RBM7

Uniprot

Q9Y580

Synonyms:

RNA-binding protein 7 (RNA-binding motif protein 7), RBM7

Immunogen:

Recombinant Human RNA-binding protein 7 protein (5-266AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A