

EMILIN1 Antibody, Biotin conjugated



PACO42709

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

May be responsible for anchoring smooth muscle cells to elastic fibers, and may be involved not only in the formation of the elastic fiber, but also in the processes that regulate vessel assembly. Has cell adhesive capacity.

Gene ID:

EMILIN1

Uniprot

Q9Y6C2

Synonyms:

EMILIN-1 (Elastin microfibril interface-located protein 1) (Elastin microfibril interfacier 1), EMILIN1, EMI

Immunogen:

Recombinant Human EMILIN-1 protein (1-347AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A