

COX11 Antibody, FITC conjugated



PACO42716

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Exerts its effect at some terminal stage of cytochrome c oxidase synthesis, probably by being involved in the insertion of the copper B into subunit I.

Gene ID:

COX11

Uniprot

Q9Y6N1

Synonyms:

Cytochrome c oxidase assembly protein COX11, mitochondrial, COX11

Immunogen:

Recombinant Human Cytochrome c oxidase assembly protein COX11, mitochondrial protein (115-276AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A