

USH1C Antibody, Biotin conjugated



PACO42721

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Required for normal development and maintenance of cochlear hair cell bundles. Anchoring/scaffolding protein that is a part of the functional network formed by USH1C, USH1G, CDH23 and MYO7A that mediates mechanotransduction in cochlear hair cells. Required for normal hearing.

Gene ID:

USH1C

Uniprot

Q9Y6N9

Synonyms:

Harmonin (Antigen NY-CO-38/NY-CO-37) (Autoimmune enteropathy-related antigen AIE-75) (Protein PDZ-73) (Renal carcinoma antigen NY-REN-3) (Usher syndrome type-1C protein), USH1C, AIE75

Immunogen:

Recombinant Human Harmonin protein (1-533AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A