

CA5B Antibody, HRP conjugated



PACO42723

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Reversible hydration of carbon dioxide.

Gene ID:

CA5B

Uniprot

Q9Y2D0

Synonyms:

Carbonic anhydrase 5B, mitochondrial (EC 4.2.1.1) (Carbonate dehydratase VB) (Carbonic anhydrase VB) (CA-VB), CA5B

Immunogen:

Recombinant Human Carbonic anhydrase 5B, mitochondrial protein (34-314AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A