

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, IHC

**Recommended dilutions:**

ELISA:1:2000-1:10000, IHC:1:20-1:200

**Protein Background:**

Component of the coenzyme Q biosynthetic pathway. May play a role in organizing a multi-subunit COQ enzyme complex required for coenzyme Q biosynthesis. Required for steady-state levels of other COQ polypeptides.

**Gene ID:**

COQ4

**Uniprot**

Q9Y3A0

**Synonyms:**

Ubiquinone biosynthesis protein COQ4 homolog, mitochondrial (Coenzyme Q biosynthesis protein 4 homolog), COQ4

**Immunogen:**


Recombinant Human Ubiquinone biosynthesis protein COQ4 homolog, mitochondrial protein (31-265AA).

**Storage:**


Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---



Immunohistochemistry of paraffin-embedded human liver tissue using PACO42726 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human spleen tissue using PACO42726 at dilution of 1:100.