

# NIP7 Antibody, Biotin conjugated



PACO42761

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

Required for proper 34S pre-rRNA processing and 60S ribosome subunit assembly.

**Gene ID:**

NIP7

**Uniprot**

Q9Y221

**Synonyms:**

60S ribosome subunit biogenesis protein NIP7 homolog (KD93), NIP7

**Immunogen:**

Recombinant Human 60S ribosome subunit biogenesis protein NIP7 homolog protein (1-180AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A