

# FARSA Antibody, HRP conjugated



PACO42763

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

FARSA

**Uniprot**

Q9Y285

**Synonyms:**

Phenylalanine--tRNA ligase alpha subunit (EC 6.1.1.20) (CML33) (Phenylalanyl-tRNA synthetase alpha subunit) (PheRS), FARSA, FARS FARSL FARSLA

**Immunogen:**

Recombinant Human Phenylalanine--tRNA ligase &alpha; subunit protein (2-508AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A