

FARSA Antibody, FITC conjugated



PACO42764

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

FARSA

Uniprot

Q9Y285

Synonyms:

Phenylalanine--tRNA ligase alpha subunit (EC 6.1.1.20) (CML33) (Phenylalanyl-tRNA synthetase alpha subunit) (PheRS), FARSA, FARS FARSL FARSLA

Immunogen:

Recombinant Human Phenylalanine--tRNA ligase α subunit protein (2-508AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A