

EFR3B Antibody



PACO42774

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

Gene ID:

EFR3B

Uniprot

Q9Y2G0

Synonyms:

Protein EFR3 homolog B, EFR3B, KIAA0953


Immunogen:

Recombinant Human Protein EFR3 homolog B protein (1-298AA).


Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human colon cancer using PACO42774 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human ovarian cancer using PACO42774 at dilution of 1:100.