

KLK5 Antibody



PACO42778

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:2000,
IHC:1:20-1:200

Protein Background:

May be involved in desquamation.

Gene ID:

KLK5

Uniprot

Q9Y337

Synonyms:

Kallikrein-5 (EC 3.4.21. -) (Kallikrein-like protein 2) (KLK-L2) (Stratum corneum tryptic enzyme), KLK5, SCTE


Immunogen:

Recombinant Human Kallikrein-5 protein (23-293AA).


Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Western blot. All lanes: KLK5 antibody at 2 μ g/ml + MCF-7 whole cell lysate. Secondary. Goat polyclonal to rabbit IgG at 1/10000 dilution. Predicted band size: 32 kDa. Observed band size: 32 kDa.



Immunohistochemistry of paraffin-embedded human skin tissue using PACO42778 at dilution of 1:100.