

KLK5 Antibody, Biotin conjugated



PACO42781

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

May be involved in desquamation.

Gene ID:

KLK5

Uniprot

Q9Y337

Synonyms:

Kallikrein-5 (EC 3.4.21. -) (Kallikrein-like protein 2) (KLK-L2) (Stratum corneum tryptic enzyme), KLK5, SCTE

Immunogen:

Recombinant Human Kallikrein-5 protein (23-293AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A