

EXOSC1 Antibody, Biotin conjugated



PACO42793

Product Information

Size:	Protein Background:
50ug	
Reactivity:	Gene ID:
Human	EXOSC1
Source:	Uniprot
Rabbit	Q9Y3B2
Isotype:	Synonyms:
IgG	Exosome complex component CSL4 (Exosome component 1), EXOSC1, CSL4
Applications:	Immunogen:
ELISA	Recombinant Human Exosome complex component CSL4 protein (9-195AA).
Recommended dilutions:	Storage:
	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A