

# TMED3 Antibody, HRP conjugated



PACO42795

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

Potential role in vesicular protein trafficking, mainly in the early secretory pathway. Contributes to the coupled localization of TMED2 and TMED10 in the cis-Golgi network.

**Gene ID:**

TMED3

**Uniprot**

Q9Y3Q3

**Synonyms:**

Transmembrane emp24 domain-containing protein 3 (Membrane protein p24B) (p24 family protein gamma-4) (p24gamma4) (p26), TMED3, C15orf22

**Immunogen:**

Recombinant Human Transmembrane emp24 domain-containing protein 3 protein (24-180AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A