

# TIMM9 Antibody, Biotin conjugated



PACO42809

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

Mitochondrial intermembrane chaperone that participates in the import and insertion of multi-pass transmembrane proteins into the mitochondrial inner membrane. May also be required for the transfer of beta-barrel precursors from the TOM complex to the sorting and assembly machinery (SAM complex) of the outer membrane. Acts as a chaperone-like protein that protects the hydrophobic precursors from aggregation and guide them through the mitochondrial intermembrane space.

**Gene ID:**

TIMM9

**Uniprot**

Q9Y5J7

**Synonyms:**

Mitochondrial import inner membrane translocase subunit Tim9, TIMM9, TIM9 TIM9A TIMM9A

**Immunogen:**

Recombinant Human Mitochondrial import inner membrane translocase subunit Tim9 protein (2-89AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A