

DHRS7 Antibody, FITC conjugated



PACO42824

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

DHRS7

Uniprot

Q9Y394

Synonyms:

Dehydrogenase/reductase SDR family member 7 (EC 1.1. -. -) (Retinal short-chain dehydrogenase/reductase 4) (retSDR4) (Short chain dehydrogenase/reductase family 34C member 1), DHRS7, DHRS7A RETSDR4 SDR34C1

Immunogen:

Recombinant Human Dehydrogenase/reductase SDR family member 7 protein (29-339AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A