

# UFC1 Antibody, Biotin conjugated



PACO42833

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

E2-like enzyme which forms an intermediate with UFM1 via a thioester linkage.

**Gene ID:**

UFC1

**Uniprot**

Q9Y3C8

**Synonyms:**

Ubiquitin-fold modifier-conjugating enzyme 1 (Ufm1-conjugating enzyme 1), UFC1

**Immunogen:**

Recombinant Human Ubiquitin-fold modifier-conjugating enzyme 1 protein (1-167AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A