

MTMR6 Antibody, Biotin conjugated



PACO42849

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Phosphatase that acts on lipids with a phosphoinositol headgroup. Acts as a negative regulator of KCNN4/KCa3.1 channel activity in CD4+ T-cells possibly by decreasing intracellular levels of phosphatidylinositol 3 phosphatase. Negatively regulates proliferation of reactivated CD4+ T-cells.

Gene ID:

MTMR6

Uniprot

Q9Y217

Synonyms:

Myotubularin-related protein 6 (EC 3.1.3. -), MTMR6

Immunogen:

Recombinant Human Myotubularin-related protein 6 protein (1-260AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A