

# FAM32A Antibody



PACO42866

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, IHC

**Recommended dilutions:**

ELISA:1:2000-1:10000, IHC:1:20-1:200

**Protein Background:**

Isoform 1, but not isoform 2 or isoform 3, may induce G2 arrest and apoptosis. May also increase cell sensitivity to apoptotic stimuli.

**Gene ID:**

FAM32A

**Uniprot**

Q9Y421

**Synonyms:**

Protein FAM32A (Ovarian tumor-associated gene 12) (OTAG-12), FAM32A, OTAG12

**Immunogen:**


Recombinant Human Protein FAM32A protein (1-112AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---



Immunohistochemistry of paraffin-embedded human ovarian cancer using PACO42866 at dilution of 1:100.