

FAM32A Antibody, Biotin conjugated



PACO42869

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Isoform 1, but not isoform 2 or isoform 3, may induce G2 arrest and apoptosis. May also increase cell sensitivity to apoptotic stimuli.

Gene ID:

FAM32A

Uniprot

Q9Y421

Synonyms:

Protein FAM32A (Ovarian tumor-associated gene 12) (OTAG-12), FAM32A, OTAG12

Immunogen:

Recombinant Human Protein FAM32A protein (1-112AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A