

TIMM8B Antibody, HRP conjugated



PACO42875

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Probable mitochondrial intermembrane chaperone that participates in the import and insertion of some multi-pass transmembrane proteins into the mitochondrial inner membrane. Also required for the transfer of beta-barrel precursors from the TOM complex to the sorting and assembly machinery (SAM complex) of the outer membrane. Acts as a chaperone-like protein that protects the hydrophobic precursors from aggregation and guide them through the mitochondrial intermembrane space.

Gene ID:

TIMM8B

Uniprot

Q9Y5J9

Synonyms:

Mitochondrial import inner membrane translocase subunit Tim8 B (DDP-like protein) (Deafness dystonia protein 2), TIMM8B, DDP2 DDPL TIM8B

Immunogen:

Recombinant Human Mitochondrial import inner membrane translocase subunit Tim8 B protein (33-83AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A