

# HOXB7 Antibody



PACO42927

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, IHC

**Recommended dilutions:**

ELISA:1:2000-1:10000, IHC:1:20-1:200

**Protein Background:**

Sequence-specific transcription factor which is part of a developmental regulatory system that provides cells with specific positional identities on the anterior-posterior axis.

**Gene ID:**

HOXB7

**Uniprot**

P09629

**Synonyms:**

Homeobox protein Hox-B7 (Homeobox protein HHO. C1) (Homeobox protein Hox-2C), HOXB7, HOX2C

**Immunogen:**


Recombinant Human Homeobox protein Hox-B7 protein (1-120AA).

**Storage:**


PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Product Images

---



Immunohistochemistry of paraffin-embedded human skin tissue using PACO42927 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human pancreatic tissue using PACO42927 at dilution of 1:100.