

# HNRNPD Antibody



PACO42991

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, IHC

**Recommended dilutions:**

ELISA:1:2000-1:10000, IHC:1:20-1:200

**Protein Background:**

**Gene ID:**

HNRNPD

**Uniprot**

Q14103

**Synonyms:**

Heterogeneous nuclear ribonucleoprotein D0 (hnRNP D0) (AU-rich element RNA-binding protein 1), HNRNPD, AUF1 HNRPD

**Immunogen:**


Recombinant Human Heterogeneous nuclear ribonucleoprotein D0 protein (1-306AA).

**Storage:**


PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Product Images

---



Immunohistochemistry of paraffin-embedded human tonsil tissue using PACO42991 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human endometrial cancer using PACO42991 at dilution of 1:100.