

# TDG Antibody



PACO43004

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

ELISA:1:2000-1:10000, WB:1:500-1:2000

**Protein Background:**

**Gene ID:**

TDG

**Uniprot**

Q13569

**Synonyms:**

G/T mismatch-specific thymine DNA glycosylase (EC 3.2.2.29) (Thymine-DNA glycosylase) (hTDG), TDG

**Immunogen:**


Recombinant Human G/T mismatch-specific thymine DNA glycosylase protein (141-410AA).

**Storage:**

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Product Images

---



Western blot. All lanes: G/T mismatch-specific thymine DNA glycosylase antibody at 8 $\mu$ g/ml + Mouse brain tissue. Secondary. Goat polyclonal to rabbit IgG at 1/10000 dilution. Predicted band size: 46 kDa. Observed band size: 46 kDa.