

IRAK3 Antibody



PACO43035

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Inhibits dissociation of IRAK1 and IRAK4 from the Toll-like receptor signaling complex by either inhibiting the phosphorylation of IRAK1 and IRAK4 or stabilizing the receptor complex.

Gene ID:

IRAK3

Uniprot

Q9Y616

Synonyms:

Interleukin-1 receptor-associated kinase 3 (IRAK-3) (EC 2.7.11.1) (IL-1 receptor-associated kinase M) (IRAK-M), IRAK3

Immunogen:

Recombinant Human Interleukin-1 receptor-associated kinase 3 protein (367-596AA).

Storage:

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images

N/A

N/A