

PACO43100

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

Binds to the consensus sequence 5'-AGAGTG-3' and has transcriptional activator activity. May have a tumor-suppressor function. May play a role in apoptosis.

Gene ID:

CSRNP1

Uniprot

Q96S65

Synonyms:

Cysteine/serine-rich nuclear protein 1 (CSRNP-1) (Axin-1 up-regulated gene 1 protein) (Protein URAX1) (TGF-beta-induced apoptosis protein 3) (TAIP-3), CSRNP1, AXUD1 TAIP3


Immunogen:

Recombinant Human Cysteine/serine-rich nuclear protein 1 protein (1-200AA).


Storage:

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Immunohistochemistry of paraffin-embedded human colon cancer using PACO43100 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human placenta tissue using PACO43100 at dilution of 1:100.