

ACTR1A Antibody



PACO43116

Product Information

Size:

50ul

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:2000

Protein Background:

Component of a multi-subunit complex involved in microtubule based vesicle motility. It is associated with the centrosome.

Gene ID:

ACTR1A

Uniprot

P61163

Synonyms:

Alpha-centractin (Centractin) (ARP1) (Actin-RPV) (Centrosome-associated actin homolog), ACTR1A, CTRN1


Immunogen:

Recombinant Human Alpha-centractin protein (187-376AA).

Storage:

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Western blot. All lanes: Alpha-centractin antibody at 7 μ g/ml + Mouse brain tissue. Secondary. Goat polyclonal to rabbit IgG at 1/10000 dilution. Predicted band size: 43 kDa. Observed band size: 43 kDa.