

FOXP2 Antibody



PACO43146

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

Transcriptional repressor that may play a role in the specification and differentiation of lung epithelium. May also play a role in developing neural, gastrointestinal and cardiovascular tissues. Can act with CTBP1 to synergistically repress transcription but CTPBP1 is not essential. Plays a role in synapse formation by regulating SRPX2 levels. Involved in neural mechanisms mediating the development of speech and language.

Gene ID:

FOXP2

Uniprot

O15409

Synonyms:

Forkhead box protein P2 (CAG repeat protein 44) (Trinucleotide repeat-containing gene 10 protein), FOXP2, CAGH44 TNRC10


Immunogen:

Recombinant Human Forkhead box protein P2 protein (416-715AA).


Storage:

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Immunohistochemistry of paraffin-embedded human prostate tissue using PACO43146 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human placenta tissue using PACO43146 at dilution of 1:100.