

PACO43155

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

Calcium-regulated non-lysosomal thiol-protease.

Gene ID:

CAPN3

Uniprot

P20807

Synonyms:

Calpain-3 (EC 3.4.22.54) (Calcium-activated neutral proteinase 3) (CANP 3) (Calpain L3) (Calpain p94) (Muscle-specific calcium-activated neutral protease 3) (New calpain 1) (nCL-1), CAPN3, CANP3, CANPL3, NCL1


Immunogen:

Recombinant Human Calpain-3 protein (582-821AA).

Storage:

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Immunohistochemistry of paraffin-embedded human colon tissue using PACO43155 at dilution of 1:100.