

DCLRE1C Antibody



PACO43162

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

DCLRE1C

Uniprot

Q96SD1

Synonyms:

Protein artemis (EC 3.1. -. -) (DNA cross-link repair 1C protein) (Protein A-SCID) (SNM1 homolog C) (hSNM1C) (SNM1-like protein), DCLRE1C, ARTEMIS ASCID SCIDA SNM1C

Immunogen:

Recombinant Human Protein artemis protein (1-300AA).

Storage:

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images

N/A

N/A