

PACO43186

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, IHC

**Recommended dilutions:**

ELISA:1:2000-1:10000, IHC:1:20-1:200

**Protein Background:**

Component of a cell-surface receptor complex that mediates cell-cell interactions between muscle precursor cells. Promotes differentiation of myogenic cells.

**Gene ID:**

BOC

**Uniprot**

Q9BWW1

**Synonyms:**

Brother of CDO (Protein BOC), BOC

**Immunogen:**


Recombinant Human Brother of CDO protein (1-125AA).

**Storage:**


PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Product Images

---



Immunohistochemistry of paraffin-embedded human small intestine tissue using PACO43186 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human kidney tissue using PACO43186 at dilution of 1:100.