

# HYAL3 Antibody



PACO43199

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, IHC

**Recommended dilutions:**

ELISA:1:2000-1:10000, IHC:1:20-1:200

**Protein Background:**

**Gene ID:**

HYAL3

**Uniprot**

O43820

**Synonyms:**

Hyaluronidase-3 (Hyal-3) (EC 3.2.1.35) (Hyaluronoglucosaminidase-3) (Lung carcinoma protein 3) (LuCa-3), HYAL3, LUCA3

**Immunogen:**


Recombinant Human Hyaluronidase-3 protein (198-417AA).

**Storage:**


PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Product Images

---



Immunohistochemistry of paraffin-embedded human breast cancer using PACO43199 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human small intestine tissue using PACO43199 at dilution of 1:100.