

ALS2 Antibody



PACO43209

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

May act as a GTPase regulator. Controls survival and growth of spinal motoneurons.

Gene ID:

ALS2

Uniprot

Q96Q42

Synonyms:

Alsin (Amyotrophic lateral sclerosis 2 chromosomal region candidate gene 6 protein) (Amyotrophic lateral sclerosis 2 protein), ALS2, ALS2CR6 KIAA1563


Immunogen:

Recombinant Human Alsin protein (1-280AA).


Storage:

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Immunohistochemistry of paraffin-embedded human breast cancer using PACO43209 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human pancreatic tissue using PACO43209 at dilution of 1:100.