

HYAL1 Antibody



PACO43225

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

May have a role in promoting tumor progression. May block the TGFB1-enhanced cell growth.

Gene ID:

HYAL1

Uniprot

Q12794

Synonyms:

Hyaluronidase-1 (Hyal-1) (EC 3.2.1.35) (Hyaluronoglucosaminidase-1) (Lung carcinoma protein 1) (LuCa-1), HYAL1, LUCA1


Immunogen:

Recombinant Human Hyaluronidase-1 protein (136-435AA).


Storage:

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Immunohistochemistry of paraffin-embedded human endometrial cancer using PACO43225 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human liver tissue using PACO43225 at dilution of 1:100.