

HYAL2 Antibody



PACO43226

Product Information

Size:

50ul

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:200-1:2000,
IHC:1:20-1:200

Protein Background:

Hydrolyzes high molecular weight hyaluronic acid, to produce an intermediate-sized product which is further hydrolyzed by sperm hyaluronidase to give small oligosaccharides. Displays very low levels of activity. Associates with and negatively regulates MST1R.

Gene ID:

HYAL2

Uniprot

Q12891

Synonyms:

Hyaluronidase-2 (Hyal-2) (EC 3.2.1.35) (Hyaluronoglucosaminidase-2) (Lung carcinoma protein 2) (LuCa-2), HYAL2, LUCA2


Immunogen:

Recombinant Human Hyaluronidase-2 protein (367-473AA).


Storage:

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Western blot. All lanes: Hyaluronidase-2 antibody at 4 μ g/ml + Mouse lung tissue. Secondary. Goat polyclonal to rabbit IgG at 1/10000 dilution. Predicted band size: 54 kDa. Observed band size: 54 kDa.



Immunohistochemistry of paraffin-embedded human kidney tissue using PACO43226 at dilution of 1:100.