

PACO43340

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

Contracts smooth muscle of the gastrointestinal and genitourinary tract, regulates growth hormone release, modulates insulin release, and may be involved in the control of adrenal secretion.

Gene ID:

GAL

Uniprot

P22466

Synonyms:

Galanin peptides [Cleaved into: Galanin; Galanin message-associated peptide (GMAP)], GAL, GAL1 GALN GLNN


Immunogen:

Recombinant Human Galanin peptides protein (20-123AA).


Storage:

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.


Product Images



Immunohistochemistry of paraffin-embedded human tonsil tissue using PACO43340 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human epityphlon tissue using PACO43340 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human adrenal gland tissue using PACO43340 at dilution of 1:100.