

BAMBI Antibody



PACO43359

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

Negatively regulates TGF-beta signaling.

Gene ID:

BAMBI

Uniprot

Q13145

Synonyms:

BMP and activin membrane-bound inhibitor homolog (Non-metastatic gene A protein) (Putative transmembrane protein NMA), BAMBI, NMA


Immunogen:

Recombinant Human BMP and activin membrane-bound inhibitor homolog protein (21-152AA).


Storage:

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.


Product Images



Immunohistochemistry of paraffin-embedded human prostate tissue using PACO43359 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human placenta tissue using PACO43359 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human colon tissue using PACO43359 at dilution of 1:100.