

STAT4 Antibody



PACO43374

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:2000,
IHC:1:20-1:200

Protein Background:

Carries out a dual function: signal transduction and activation of transcription. Involved in IL12 signaling.

Gene ID:

STAT4

Uniprot

Q14765

Synonyms:

Signal transducer and activator of transcription 4, STAT4


Immunogen:

Recombinant Human Signal transducer and activator of transcription 4 protein (569-748AA).


Storage:

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Western blot. All lanes: Signal transducer and activator of transcription 4 antibody at 7 μ g/ml + Jurkat whole cell lysate. Secondary. Goat polyclonal to rabbit IgG at 1/10000 dilution. Predicted band size: 86 kDa. Observed band size: 86 kDa.



Immunohistochemistry of paraffin-embedded human cervical cancer using PACO43374 at dilution of 1:100.