

SERPINB9 Antibody



PACO43377

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:2000

Protein Background:

Granzyme B inhibitor.

Gene ID:

SERPINB9

Uniprot

P50453

Synonyms:

Serpin B9 (Cytoplasmic antiproteinase 3) (CAP-3) (CAP3) (Peptidase inhibitor 9) (PI-9), SERPINB9, PI9


Immunogen:

Recombinant Human Serpin B9 protein (1-376AA).

Storage:

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Western blot. All lanes: Serpin B9 antibody at 3 μ g/ml + K562 whole cell lysate. Secondary. Goat polyclonal to rabbit IgG at 1/10000 dilution. Predicted band size: 42 kDa. Observed band size: 42 kDa.