

# UTS2 Antibody



PACO43382

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC

**Recommended dilutions:**

ELISA:1:2000-1:10000, WB:1:200-1:1000,  
IHC:1:20-1:200

**Protein Background:**

Highly potent vasoconstrictor.

**Gene ID:**

UTS2

**Uniprot**

O95399

**Synonyms:**

Urotensin-2 (Urotensin II) (U-II) (UII), UTS2

**Immunogen:**


Recombinant Human Urotensin-2 protein (35-124AA).

**Storage:**


PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Product Images

---



Western blot. All lanes: Urotensin-2 antibody at 8 $\mu$ g/ml + Rat brain tissue. Secondary. Goat polyclonal to rabbit IgG at 1/10000 dilution. Predicted band size: 15, 17 kDa. Observed band size: 15 kDa.



Immunohistochemistry of paraffin-embedded human pancreatic tissue using PACO43382 at dilution of 1:100.