

PRKG1 Antibody



PACO43397

Product Information

Size:

50ul

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:2000

Protein Background:**Gene ID:**

PRKG1

Uniprot

Q13976

Synonyms:

cGMP-dependent protein kinase 1 (cGK 1) (cGK1) (EC 2.7.11.12) (cGMP-dependent protein kinase I) (cGKI), PRKG1, PRKG1B PRKGR1A PRKGR1B


Immunogen:

Recombinant Human cGMP-dependent protein kinase 1 protein (1-300AA).

Storage:

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Western blot. All lanes: cGMP-dependent protein kinase 1 antibody at 3 μ g/ml + Mouse lung tissue. Secondary. Goat polyclonal to rabbit IgG at 1/10000 dilution. Predicted band size: 77, 78, 45 kDa. Observed band size: 77 kDa.