

PACO43407

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

This enzyme is specific for acyl chain lengths of 4 to 16.

Gene ID:

ACADM

Uniprot

P11310

Synonyms:

Medium-chain specific acyl-CoA dehydrogenase, mitochondrial (MCAD) (EC 1.3.8.7),
ACADM


Immunogen:

Recombinant Human Medium-chain specific acyl-CoA dehydrogenase, mitochondrial
protein (26-421AA).

Storage:

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Immunohistochemistry of paraffin-embedded human kidney tissue using PACO43407 at dilution of 1:100.